



AEROBIOTIX™

BETTER AIR FOR HEALTHCARE

Aerobiotix Illuvia 500 Specification sheet

Physical Dimensions

- 18 x 18 x 54 inches
- Weighs 150 pounds

Electrical Ratings

- Rated at 110 volts
- Rated at 10 Amps
- Grounded outlet required

Decibel Level

- 62-64 decibels

Certifications

- UL 507
- 60601-1-2 EMC

Air Circulation

- 450 CFM
- Will not affect positive pressure
 - Air intake matches air exhaust
- Will not affect room air exchanges
 - In room air cleaners approved by ASHRE and CDC for use in OR
- Operating Room Flow
 - Will not affect room flow, exhaust air is baffled, and vented over a larger surface area.

Cleaning

- Wipe down entire machine with antibiotic wipe.

